

Name	Title	Date	Unit
Stryker S3 bed – Skills Checklist			✓
1. Siderails and Brakes <ul style="list-style-type: none"> • Siderails up-down * <i>Use intermediate position with patients getting in/out of bed</i> • Brake 			
2. UP/DOWN Functions <ul style="list-style-type: none"> • Bed up/down • Fowlers up/down • Knee up/down • Cardiac Chair 			
3. Foot End Controls <ul style="list-style-type: none"> • Lockout functions • Trendelenburg/Reverse 			
4. Scale <ul style="list-style-type: none"> • Zero • Weigh 			
5. Menu functions <ul style="list-style-type: none"> • lbs/kg • alarm tone 			
6. Chaperone Bed Exit System <i>To Arm System</i> <ol style="list-style-type: none"> 1. Zero the Scale without patient in the bed 2. Arm – Push and release the Arm/Disarm key 3. Select Zone- push Zone key until desired zone lights 4. To Disarm- push and hold Disarm key 			
7. iBED Awareness <i>Set bed configuration</i> <ol style="list-style-type: none"> 1. Push and release iBED key 2. All 3 green lights on footboard If changed- amber light will flash, - display and dashboard will identify change 3. To Disarm iBED Push and hold iBED key lights will turn off 			
8. Describe scenario where iBed awareness function would apply			
9. Manual CPR			
10. Steering			
11. Foley bag hooks			
12. Cords- plug in electrical and Nurse Call cable			
13. Identify electrical outlet location on foot end of the bed			
Skills Validator (Print)		Skill Validator (Signature)	
Employee Name (Print)		Employee Signature	